

[³⁵S]-amino acids

³⁵ S	L-Methionine, [³⁵S]-	CH ₃ S(CH ₂) ₂ CH(NH ₂)COOH; In 50 mM Tris-HCl, pH 7.4 containing 10 mM 2-mercaptoethanol; > 1000 Ci/mmol; 10-20 mCi/ml
	Met-^[35S]-label	L-Methionine, [³⁵ S], L-Cysteine, [³⁵ S]; E. coli hydrolysate containing 70% L- Methionine, [³⁵ S]-, ~25% L-Cysteine, [³⁵ S]- and ~5% miscellaneous amino acids; > 1000 Ci/mmol
	L-Cystein, [³⁵S]-	HSCH ₂ CH(NH ₂)COOH; In stabilized aqueous solution containing 10 mM DTT; >1000 Ci/mmol; 10 mCi/ml

[³⁵S]-nucleotide analogue

³⁵ S	3'-Phosphoadenosine-5'- phosphosulfate (PAPS), [³⁵S]-	In aqueous solution containing 50% ethanol; >1000 Ci/mmol; 1 mCi/ml
-----------------	---	---

You need different concentrations or specific activities. Please send us your request.